

# MORGANVALE 2020

| LOT NO | ANIMAL IDENT | WT   | AGE       | SS | E/PIG   | POLL GENE | HEIFER BULLS |
|--------|--------------|------|-----------|----|---------|-----------|--------------|
| 1      | WLMP045      | 1030 | 22.5      | 42 | 100/100 | PP        |              |
| 2      | WLMP104      | 888  | 22        | 43 | 100/90  | PP        |              |
| 3      | WLMP017      | 1050 | 23        | 41 | 100/100 | PP        |              |
| 4      | WLMP071      |      | Withdrawn |    |         |           |              |
| 5      | WLMP027      | 1000 | 23        | 40 | 100/100 | PP        |              |
| 6      | WLMP082      | 1020 | 22        | 41 | 100/100 | PP        |              |
| 7      | WLMP059      | 916  | 22        | 43 | 100/100 | PP        |              |
| 8      | WLMP014      | 944  | 23        | 45 | 100/100 | PP        | H            |
| 9      | WLMP073      | 984  | 22        | 43 | 100/100 | PP        | H            |
| 10     | WLMP047      | 960  | 22.5      | 44 | 100/100 | PP        |              |
| 11     | WLMP008      | 918  | 23        | 43 | 100/100 | PP        | H            |
| 12     | WLMP042      | 894  | 23        | 40 | 100/100 | PP        | H            |
| 13     | WLMP016      | 918  | 23        | 43 | 100/100 | PP        | H            |
| 14     | WLMP015      | 924  | 23        | 40 | 100/100 | PP        |              |
| 15     | WLMP174      | 874  | 20        | 41 | 100/100 | PP        |              |
| 16     | WLMP078      | 894  | 22        | 40 | 100/100 | PP        |              |
| 17     | WLMP151      | 880  | 21        | 40 | 100/100 | PP        |              |
| 18     | WLMP055      | 908  | 22.5      | 42 | 100/100 | PP        | H            |
| 19     | WLMP235      | 892  | 19        | 38 | 100/100 | PP        |              |
| 20     | WLMP246      | 882  | 19        | 44 | 100/100 | PP        |              |
| 21     | WLMP030      |      | Withdrawn |    |         |           |              |
| 22     | WLMP161      |      | Withdrawn |    |         |           |              |
| 23     | WLMP158      | 904  | 20.5      | 41 | 100/100 | PP        |              |
| 24     | WLMP183      | 910  | 20        | 42 | 100/100 | PP        |              |
| 25     | WLMP219      | 886  | 19.5      | 38 | 100/100 | Ph        |              |
| 26     | WLMP218      | 888  | 19.5      | 40 | 100/100 | PP        |              |
| 27     | WLMP169      | 904  | 20        | 42 | 100/100 | PP        |              |
| 28     | WLMP129      | 886  | 21        | 37 | 100/100 | PP        |              |
| 29     | WLMP260      | 772  | 19        | 38 | 90/90   | PP        |              |
| 30     | WLMP005      | 868  | 23        | 43 | 100/100 | PP        | H            |

|    |          |     |      |    |         |            |   |
|----|----------|-----|------|----|---------|------------|---|
| 31 | WLMP236  | 856 | 19   | 41 | 100/90  | PP         |   |
| 32 | WLMP152  | 894 | 21   | 43 | 100/100 | PP         |   |
| 33 | WLMP105  | 900 | 22   | 44 | 100/100 | PP         |   |
| 34 | WLMP011  | 862 | 23   | 43 | 100/90  | PP         | H |
| 35 | WLMP136  | 848 | 21   | 40 | 100/90  | PP         | H |
| 36 | WLMP196  | 886 | 20   | 40 | 95/95   | PP         |   |
| 37 | WLMP241  | 834 | 19   | 41 | 100/100 | PP         |   |
| 38 | WLMP224  | 840 | 19.5 | 44 | 100/100 | PP         |   |
| 39 | WLMP189  | 840 | 20   | 38 | 100/100 | PP         | H |
| 40 | WLMP292  | 840 | 18.5 | 42 | 100/100 | PP         |   |
| 41 | WLMP272  | 824 | 19   | 40 | 100/100 | PP         |   |
| 42 | WLMP269  | 844 | 19   | 38 | 100/100 | PP         |   |
| 43 | WLMP166  | 822 | 20.5 | 41 | 100/100 | PP         |   |
| 44 | WLMP214  | 884 | 19.5 | 41 | 100/100 | PP         |   |
| 45 | WLMP311  | 704 | 18   | 38 | 100/100 | PP         |   |
| 46 | WLMP157  | 788 | 20.5 | 40 | 100/100 | PP         | H |
| 47 | WLMP264  | 768 | 19   | 40 | 100/100 | PP         |   |
| 48 | WLMP345  | 792 | 18.5 | 38 | 100/100 | PP         |   |
| 50 | MGVPP039 | 864 | 23.5 | 42 |         | NOT TESTED |   |
| 51 | MGVPP182 | 794 | 20   | 39 |         | NOT TESTED |   |
| 52 | MGVPP195 | 748 | 19.5 | 39 |         | NOT TESTED |   |
| 53 | MGVPP212 | 746 | 19.5 | 39 |         | NOT TESTED |   |
| 54 | MGVPP304 | 754 | 18   | 41 |         | NOT TESTED |   |

Bulls for Heifer Joinings Lots- 8 , 9 , 11, 12, 13, 18, 34, 35

Bulls for Heifer Joinings in well grown Heifers Lots- 30, 39, 46

POLL TEST PP- HOMOZYGOUS

Ph HETEROZYGOUS

All bulls weighed 19-2-20  
Drenched Dectomax Injectable 3/02/2020